

Document Detail

	RCI IBTS CM IBTS/RCI/CM/0001[4]
Title:	CUSTOMER MANUAL FOR THE RED CELL
	IMMUNOHAEMATOLOGY LABORATORY
Owner:	QA DOC CON QA DOC CONTROL
Status	CURRENT
Effective Date:	24-Feb-2025
Expiration Date:	24-Feb-2027

<u>Review</u>

Review: IBTS DOC REVIEW AND APPROVAL

Level	Owner Role	Actor	Sign-off By
1	DOCUMENT CONTROLLER	REBECCA WALDEN	REBECCA WALDEN
2	RCI WRITER IBTS	EDEL SCALLY	EDEL SCALLY
3	RCI REVIEWER IBTS	CONOR GALVIN	CONOR GALVIN
4	MEDICAL HEAD OF DEPT NBC	KIERAN MORRIS	KIERAN MORRIS
4	RCI HEAD OF DEPT NBC	EDEL SCALLY	EDEL SCALLY
4	LABS HEAD OF TESTING IBTS	PADRAIG WILLIAMS	PADRAIG WILLIAMS
5	Labs Patient Testing Quality & Regulatory Reviewer IB	MARIA KANE	
I	Review: IBTS DOC REVIEW AND APPROVAL		
<u>Level</u>	Owner Role	Actor	<u>Sign-off By</u>
Level	Owner Role DOCUMENT CONTROLLER	Actor BECKY WHITE	<u>Sign-off By</u> BECKY WHITE
<u>Level</u> 1 2			
1	DOCUMENT CONTROLLER	BECKY WHITE	BECKY WHITE
1 2	DOCUMENT CONTROLLER RCI WRITER IBTS	BECKY WHITE EDEL SCALLY	BECKY WHITE EDEL SCALLY
1 2 3	DOCUMENT CONTROLLER RCI WRITER IBTS RCI REVIEWER IBTS	BECKY WHITE EDEL SCALLY LAYLA SPARLING	BECKY WHITE EDEL SCALLY LAYLA SPARLING
1 2 3 4	DOCUMENT CONTROLLER RCI WRITER IBTS RCI REVIEWER IBTS LABS HEAD OF TESTING IBTS	BECKY WHITE EDEL SCALLY LAYLA SPARLING PADRAIG WILLIAMS	BECKY WHITE EDEL SCALLY LAYLA SPARLING PADRAIG WILLIAMS

Change Orders

Changes as described on Change Order:

Change Order No.

Change Orders - Incorporated

Changes as described on Change Order:

Change Order No. IBTS/CO/0048/25

IBTS/RCI/CM/0001	Ver. 4	Page 2 of 51

TITLE: CUSTOMER MANUAL FOR THE RED CELL IMMUNOHAEMATOLOGY LABORATORY

Change Description:

- 1. Remove reference to the RCI laboratory being INAB accredited
- 2. Add participation of RCI in WASP EQA scheme
- 3. Change name of Lab Quality to Aurevia
- 4. In section 10.2 refer to IBTS/QA/SOP/0063 Attachment 6.2 on the giveblood.ie website
- 5. In Section 11 include in 'Suitable for transfusion until' section of the compatibility label is for internal use only.
- 6. Detail the cost list of services is emailed to each hospital at the start of each year or on request.
- 7. Change training type from Read & Understand to Read Only
- 8. Update Section 7 for clarity on referral test services & add anti-G investigation
- 9. Add information in compliments for patient, user and personnel feedback form
- 10. Add detail on verbal request procedure.
- 11. Add detail that samples should be received by day 4 phlebotomy to ensure there is sufficient time to process the sample
- 12. Add additional information for the processing of urgent samples.
- 13. Information added on sample rejection monitoring.
- 14. To DAT test result add detail that there is no commercial IQC available for IgA and IgM portion of the monospecific DAT
- 15. Update sections where IBTS/RCI/LM/0001 is referred out to
- 16. Remove Urine Inhibition Test

Reason for Change:

- 1. As per CC 050/25
- 2. As per CC 071/22
- 3. As per CC 011/25
- 4. For ease of access for users
- 5. Clarify the expiry time of the crossmatch is documented on the report not on the label
- 6. To provide additional information for users
- 7. As per IR 299/24 INAB reference PD/2024-01/04
- 8. To make it clear to users that samples may be referred elsewhere when the test is not performed in RCI or further testing is required elsewhere to support findings. Anti-G test provided.
- 9. Requirement for ISO 15189:22 standard Section 8.6.2 CC 308/24
- 10. Requirement of ISO 15189:22 standard Section 7.3.2.2 CC 308/24
- 11. Provide user with information on cut off times for samples
- 12. Requirement of ISO 15189:22 standard Section 7.2.6.1 CC 308/24
- 13. Requirement of ISO 15189:22 standard Section 7.2.7.3 CC 308/24
- 14. Make users aware there is no control for this portion of the DAT result. IR 315/24 INAB ref: SN-2024-1-05
- 15. New format of IBTS/RCI/LM/0001 written to adhere to ISO 15189:22 standard
- 16. As per CC 0197/24

IBTS/RCI/CM/0001	Ver. 4	Page 3 of 51
------------------	--------	--------------

Change order No.:

IBTS/CO/0048/25

Referenced Documents

IBTS/MED/SOP/0050	IBTS/RCI/LM/0001	BT - 0345
IBTS/QA/QM/0001	IBTS/RCI/SOP/0013	IBTS/RCI/FORM/0001
IBTS/EXT/DOC/0025	IBTS/EXT/DOC/0012	[IBTS/EXT/DOC/0017
IBTS/EXT/DOC/0033		20.

SmartSolve Roles

		and i
SmartSolve Roles		, S ^{io}
MED CON BMR HLA NBC	RCI THOD NBC	RCI SMS NBC
MED CON IH NBC	RCI MS NBC	RCI TMS NBC
MED CON MSD IBTS	RCI LA NBC	MED SMO NBC
Training Type		
Chaff has in adding a new income		

Training Type

Staff trained in previous versions	New staff & staff returning from long term
	leave
Read only	Read only
1P2	
SmartSolve Document Category	

SmartSolve Document Category

Category	Mobile	Sector Cryobiology	Website	GDP
Yes / No	No	No	Yes	No
	5	1	ł	
	Se.			
	23			
	11.			
Jerift	er			
2				
$P_{B_{i}}$				
1er				
~				

IBTS/RCI/CM/0001	
------------------	--

TABLE OF CONTENTS

1	INTR	ODUCTION	6	
2	QUA	LITY POLICY	8	
3	GENE 3.1 3.2	ERAL INFORMATION RCI Laboratory Laboratory Director	>	10
	3.3	Service Operating Times		10
	3.4 3.5	Key Personnel and Contact Details Sample Testing Schedule		11
	3.6	Turnaround Times		11 14
	3.7	Turnaround Times Service Fees and Charges		14
	3.8	Data Protection		14
4	REOI	JEST FORMS	14	
т	4.1	General Information : Samples and Forms	······ ± +	14
	4.2	RCI Laboratory Request Forms		16
	4.3	Ordering IBTS Request Forms		16
	4.4	Completion of Request Forms		17
5	ςδη	PLE INFORMATION	18	
5	5.1	Sample Collection		18
	5.2	Service Requestor Responsibilities		
	5.3	Sample Labelling		
6	SVV	PLE DELIVERY, PACKAGING AND TRANSPORT	23	
0	6.1	Sample Delivery		23
	6.2	Sample Packaging and Transport		
7	TEST			
-	7.1	Services Provided		26
	7.2	Referral Test Services		
	7.3	Provision of Rare Donor Red Cell Components from International Rare Blood Progra	ammes	43
	7.4	Concessionary Release of Blood Components		43
	7.5	Medical and Scientific Consultancy Service		
	7.6	Haemovigilance Advisory Service		
	7.7	Repeat Examination		
	7.8	Further Examination of the Primary Sample		14
8		RNAL AND INTERNAL QUALITY ASSURANCE SCHEMES		
	8.1	External Quality Assessment Programmes (EQA)		
	8.2	Inter-Laboratory Comparison Scheme		
	8.3 8.4	Internal Quality Assessment Programme		
	8.4 8.5	Authorisation Non-Conformance / Failure in an EQA Scheme		
0				13
9		ORTING OF RESULTS		
	9.1 9.2	Approval of Test Results and Issuing Reports		
	9.2 9.3	Issuing Reports on Critical Samples where Results are Delayed Reporting of Results by Fax / Encrypted Email		
	9.5 9.4	Telephoned Results		
	9.5	Archiving of Patients Records		
10		OMER SERVICE / SATISFACTION AND REVIEW		
10	10.1	Service Level Agreements		18
	10.1	Customer Complaints / Compliments		
		r		-

IBTS/RCI/CM/0001	Ver. 4	Page 5 of 51	
------------------	--------	--------------	--

	10.3 10.4 10.5 10.6	Quality Management Review Customer Liaison Hospital Transfusion Committees Continuous Improvement	48 49
	/ERSE EV 11.1 11.2	BILITY AND REPORTING OF SERIOUS ADVERSE REACTIONS (SARS) AND SERIOUS ENTS (SAES)OTHER REQURIREMENT	49 50
12	REFERE	NCES	
13	ATTACH	IMENTS	
	Jerift	NCES	

IBTS/RCI/CM/0001	Ver. 4	Page 6 of 51

TITLE: CUSTOMER MANUAL FOR THE RED CELL IMMUNOHAEMATOLOGY LABORATORY

1 INTRODUCTION

- **1.1** This manual is designed to provide an overview of the services available from the Red Cell Immunohaematology (RCI) Laboratory at the National Blood Centre. It is intended for the customers of both the routine compatibility services, routine antenatal services and of the referral immunohaematology service.
- **1.2** The services described are provided to hospitals, hospital blood transfusion laboratories and medical practitioners in the public and private health care sectors in the Republic of Ireland.
- **1.3** This manual specifies the minimum requirements for the labelling of samples and for the completion of request forms to ensure sufficient information is received for the requested service to be optimally delivered. The RCI laboratory will accept correctly completed request forms from designated facilities requesting its service provided the pertinent details are completed and the samples accompanying the forms meet the current specified criteria. Authorised personnel will review test request documentation to determine suitability of tests requested.
- 1.4 The RCI Laboratory is subject to regular scheduled inspection by the Health Products Regulatory Authority (HPRA) for compliance as a Blood Establishment to the relevant EU Directives and Irish Statutory Instruments (SI 360 of 2005, SI 547 of 2006, SI 562 of 2006).

1.5 The laboratory complies with S1 547 of 2006 incorporating Articles 14 and 15 of Directive 98/ 2002/EC (Traceability Requirements, Notification of SAR/E).

- **1.6** All work is carried out within the framework of a documented quality system, according to Good Laboratory Practice (GLP) and Good Manufacturing Practice (GMP). The laboratory operates to internal policies and procedures for all activities as defined by the IBTS Quality Management System. This manual is a controlled document as part of that System. All red cell referral & compatibility services undergo continuous review through quality assurance and audit activities.
- **1.7** Samples are disposed of by the laboratory in accordance with IBTS Health and Safety procedures and, in compliance with waste management regulations.

IBTS/RCI/CM/0001	Ver. 4	Page 7 of 51

- **1.8** This manual should be read in conjunction with the IBTS product master files and the RCI laboratory manual. The laboratory manual and product master files may be viewed on the giveblood.ie website.
- **1.9** IBTS laboratory management is committed to the provision of a full and effective service. To this end it ensures:
 - Optimum staff recruitment, training, development and retention at all levels.
 - Procurement, validation and maintenance of appropriate equipment /resources.
 - Maintaining sample integrity and thereby the correct performance of laboratory examinations.
 - Use of examination procedures that are fit for purpose and ensure the highest achievable quality.
 - Timely, confidential, accurate and clinically useful reporting of examination results.
 - Assessment of customer satisfaction, in addition to internal audit and external quality assessment.
 - Notification to users of significant changes to RCI laboratory processes/procedures where the results or their interpretation could be significantly different, prior to implementation.

1.10 A copy of this manual is available on the internet at: <u>https://www.giveblood.ie/Clinical-Services/Red-Cell-Immunohaematology-Diagnostics/</u>. Hard copies of the customer guide will not be supplied.

- **1.11** When key changes are made to either the tests or the services identified in this manual, the customer will be notified in writing or by email. The electronic copy of the manual will be modified and made available to the customer on the website above.
- **1.12** The term 'BSH Guidelines 2012' shall refer to 'Guidelines for pre-transfusion compatibility Procedures in Blood Transfusion Laboratories' British Committee for Standards in Haematology, 2012, throughout the document.
- **1.13** The term 'Hospital Blood Transfusion Laboratory' is used to describe the Blood Transfusion Laboratories, in hospitals to which the RCI laboratory provides a referral service.
- **1.14** The term 'Hospital Blood Bank' is used to describe the situation where the RCI laboratory act as an institution's Hospital Blood Bank.

2 QUALITY POLICY

The RCI management team will ensure the laboratory complies with the IBTS Quality Policy laid out in *IBTS/QA/QM/0001*. In addition the RCI Laboratory will comply with the standards set by, AML-BB, S.I. No. 360 of 2005 European Communities (Quality and Safety of Human Blood and Blood Components) Regulations 2005 and EU Directive 2002/98/EC (setting standards of quality and safety for the collection, testing, processing, storage and distribution of human blood and blood components) for the services and tests defined in this manual and is committed to:


- Staff recruitment, training and development at all levels to provide an effective and efficient service to its users.
- Providing and managing resources to ensure that all examinations are processed to produce the highest quality results possible.
- Reporting results in ways, which are timely, confidential, accurate and are supported by clinical advice and interpretation when required,
- Implementation of Internal Quality Control, External Quality Assessment, Audit and Assessment of Customer Satisfaction to continuously improve the quality of the service
- Compliance with relevant environmental legislation.
- Adherence to appropriate technical and professional standards.
- Management and staff are committed to creating a quality culture within the Department by continuously improving our services based on the results of performance through data review, internal quality audits, equipment maintenance, Quality Control programmes and the assessment of customers' needs.

IBTS/RCI/CM/0001	Ver. 4	Page 9 of 51

3 GENERAL INFORMATION

3.1 RCI Laboratory

Location: The RCI Laboratory is based at the National Blood Centre (NBC) in Dublin. The NBC is located on the site of St James Hospital (see map below):



IBTS/RCI/CM/0001	Ver. 4	Page 10 of 51
------------------	--------	---------------



3.2 Laboratory Director

The RCI Laboratory is directed by Consultant Haematologist; Dr Kieran Morris.

3.3 Service Operating Times

Department / Activity	Opening Hours RCI
Routine Laboratory Emergency Out of Hours Service* (On-Call Service)	Monday to Friday 8.00 am – 19:00 pm Excluding Bank Holidays Monday to Friday: 19:00 pm to 8:30 am (Scientist On-call**) Saturday, Sundays & Bank Holidays:
Sample Reception	Scientist On-Call 24 hrs.** Security: 24 Hrs

- * The out-of-hours service is for emergency referrals only and where immediate blood component transfusion is clinically indicated.
- ** The Scientist on-call for the RCI Laboratory is not on-site & is called in as/when required.

IBTS/RCI/CM/0001	Ver. 4	Page 11 of 51	
------------------	--------	---------------	--

3.4 Key Personnel and Contact Details

SECTION	RCI, NBC
Consultant Haematologist	Dr Kieran Morris 01-4322800
	or Specialist Medical Officer / Registrar on duty 01-4322800
Chief Medical Scientist	Ms Edel Scally 01-4322966
Laboratory (Routine Hours)	01-4322972 01-4322973
Laboratory (Out of Hours)	01-4322800 (Switch) Ask for Medical Scientist On-call.
Clinical issues (Out of Hours)	01-4322800 (Switch) Ask for doctor on duty/call.
Components Issue Laboratory <u>Hospital Blood Bank Service</u> <u>ONLY</u> *	Routine Hours (08:30 am – 19:00 pm): Contact Hospital Services (01-4322970) (Hospital Services staff will contact staff in the RCI Laboratory to fill the order)
Laboratory Fax No.	Out-of-hours: 01 4322800 (Switch) 01-4322709
Switch	01-4322800
Emergency Contact No. (Hospital Services Department)	01 4540131

*All other platelet orders (i.e. from Hospital Blood Transfusion Laboratories) are handled by the Components Issue Laboratory

3.5 Sample Testing Schedule

3.5.1 Routine Service

3.5.1.1 Hospital Blood Bank Service

The RCI laboratory will accept and process samples for routine compatibility testing throughout the day Monday to Friday. Sample testing on any day is prioritised based on clinical need. The cut off time for receipt of samples for group & hold / routine crossmatch on the same day is 5pm.

3.5.1.2 Hospital Blood Transfusion Laboratories (referral service)

Samples should be received by RCI by day 4 post bleed to allow for sufficient time to perform the required testing. If the sample is >4 days old on receipt, the laboratory will determine whether the sample

can be accepted or not. The referring laboratory will be contacted if the sample is rejected.

Samples are processed based on clinical need (except for batched tests – see below). Samples which are referred for antibody investigation and provision of blood will be prioritised for testing during the routine working day. The cut off time for sample receipt for provision of blood during the routine working day is 2pm (please ensure that samples for crossmatching are sent without delay and directly to the laboratory to meet this cut off time).

In an urgent situation the RCI laboratory should be contacted by telephone and provision will be made to process the sample urgently or out of hours if required. Samples for antenatal antibody titration/ quantitation will be batched and processed to meet the test turnaround times specified in Section 7.1.1.

3.5.2 Emergency Service

3.5.2.1 Hospital Blood Bank Service

The RCI laboratory provides an emergency blood group and compatibility service, both routine and out of hours, for organisations where the IBTS acts as their Hospital Blood Bank and have in place a service level agreement (SLA). The RCI Laboratory maintains a supply of 2 units of O RhD negative, un-crossmatched blood at the NBC for issue in emergency issue to the Royal Victoria Eye & Ear Hospital. In this instance the units are labelled as Emergency Stock on the BT -0396 Compatibility label.

3.5.2.2 Referral Service for Hospital Transfusion Laboratories

The RCI Laboratory also provides an emergency immunohaematology / compatibility testing referral service for Hospital Blood Transfusion Laboratories.

3.5.2.3 Requesting Emergency Services

If a sample is urgent please indicate this on the request form by completing the 'Treat as Emergency' box to ensure that the request is prioritised by the laboratory.

Please contact the laboratory to discuss the urgent requirement.

• Routine Hours: Contact the laboratory directly (refer to Section 3.4 for contact details).

Out of hours: Contact switch and request to speak to the medical personnel on-call (refer to Section 3.4 for contact details).
 Requests for emergency services / compatibility testing out of hours will be assessed in accordance with the urgency of the request by the IBTS Specialist Medical Officer (SpMO) / Registrar / Consultant Haematologist.

3.5.2.4 Procedure for Urgent Requests

All urgent requests must be phoned to the laboratory and the priority of the sample indicated on the BT-345 request form. Urgent samples must be transported as urgent with clear instructions on delivery of the sample. The turnaround time for the sample will be as indicated in Section 7 below, however on receipt of the sample will be prioritised. Where multiple urgent samples are received, samples will be triaged by the Medical Team.

When making urgent request the following details will be requested and confirmed (by reading back to the person giving the information):

- Hospital / ward
- Name of person making the request and contact details
- The urgency of the request (date and time required) and estimated time of sample arrival
- Patient's name (if known), hospital / emergency / trauma number and date of birth.

In addition the following details will be confirmed to the medical scientist:

- Number and type of blood component requested
- Blood Group, ABO/ Rh/ K type, if known (from referral laboratory only), serology results at the referring site & details of known antibodies
- Reason for transfusion
- Transfusion history (if known)
- Relevant clinical condition
- Current haemoglobin

Note:

1	. Where a clinical condition dictates that a transfusion is required prior to the completion of testing, the transfusion support may vary depending on the degree of clinical urgency, the availability of an emergency stock of red cells on site at the hospital and prior availability of the patient's sample and validated
2	 blood group at the RCI laboratory. Where blood is required urgently but prior to the completion of compatibility testing and the patient's sample is not known to contain clinically significant antibodies, transfusion support will be provided by the RCI laboratory or issued
3	 from the hospital's own stock. Where antibodies are detected in the patient's sample, the relative risks of abbreviated testing prior to emergency transfusion will need to be discussed between the IBTS Consultant Haematologist/ Registrar/ SpMO and the clinician
4	 in charge at the hospital. Where there is concessionary release of a product/component or a deviation from standard procedure a concessionary release will be authorised by the IBTS

IBTS/RCI/CM/0001	Ver. 4	Page 14 of 51
IBTS/RCI/CM/0001	Ver. 4	Page 14 of 51

Consultant Haematologist/ Registrar/SpMO following consultation with the patient's attending clinician (in accordance with IBTS/MED/SOP/0050).
If a verbal request is received for additional units or testing the referring hospital must send a completed BT-345 or BT – 0007 request form for the additional testing required. Where a verbal request for testing is received, testing may commence prior to receipt of the request form but a request form must be received for every requested test.

3.6 Turnaround Times

For turnaround times for routine and emergency requests see Section 7.1.1.

3.6.1 Non-compliance with turnaround times

- Turnaround times are monitored monthly.
- Should there be a significant delay in the expected turnaround times, the requestor will be notified in the instances where the delay could compromise patient care.
- The requesting facility must inform the laboratory of any change in the urgency of the blood so that appropriate action can be taken.

Note:

Overuse of the urgent service will adversely affect the turnaround time of all urgent tests.

3.7 Service Fees and Charges

A list of current services fees and charges, with detailed information regarding out of hours service charges, are available from: Director of Finance - 01 4322800. The cost list of service is emailed to all hospitals at the beginning of each year and is available from the IBTS on request.

3.8 Data Protection

Under the General Data Protection Regulation (GDPR) (EU) 679/16 and the Data Protection Act 2018, the RCI Laboratory acts as a DATA PROCESSOR for the referring laboratory/organisation when samples are referred for testing to IBTS. This testing provides diagnostic testing for patients of the referring laboratory. The referring laboratories are the DATA CONTROLLERS.

4 REQUEST FORMS

4.1 General Information : Samples and Forms

- It is the policy of the IBTS laboratories to treat all samples as potentially infectious or high risk. Therefore, we advise that universal precautions be taken in the collection, packaging and the delivery of samples being sent to the laboratories for analysis.
- All materials used in the collection of samples should be treated as

potentially hazardous and discarded according to the hospital guidelines for waste management and in compliance with relevant regulations.

- Samples for referral should be <u>freshly drawn</u> venous samples without dilution by intravenous fluid. Referred samples should not have been tested/sub-sampled at the referring hospital; exceptions can be made for patients that are difficult to sample e.g. poor veins, following discussion with the RCI laboratory.
- All patient samples must be labelled at the bedside applying positive patient identification.
- Sending haemolysed and/or lipaemic samples should be avoided where possible as free haemoglobin and/or fatty plasma can produce test result errors (especially when using automated equipment). Such samples may have to be rejected. However, it is recognised that there are situations when haemolysis, in particular, is a result of the patient's condition.
- Samples referred to the RCI laboratory should conform to the requirements for the timing of sample collection, as defined in Section 5.3.6.
- Samples should be transported promptly to the RCI laboratory (Samples should not be stored overnight in transport vehicles).
- Samples must not be exposed to direct sunlight or extremes of temperature; samples should be transported in an ambient temperature range (2 – 25°C) unless otherwise specified in Section 7.1.1. If samples are required to be stored prior to referral to the RCI laboratory, they should be stored at 2-8°C (exceptions are details in Section 7.1.1). This is the responsibility of the referring hospital / institution.
- Sample forms / packaging are date and time stamped on receipt at the IBTS.
- On receipt in the laboratory, samples are registered with a unique RCI sample number and all stored aliquots from the primary sample are labelled with the assigned RCI sample number.
- RCI scientific staff will review request forms and samples against predefined acceptance criteria to determine if they are suitable for the tests requested (See Sections 4.4 and 5.3 for request form and specimen labelling requirements, respectively) Where it is determined that the request form and/or sample is not suitable, the requesting hospital will be informed.

Note:

1.	Incorrect or incomplete forms/ samples may result in the tests not
	being undertaken and may require a second sample to be submitted
	resulting in increased turnaround times with service delay.
2.	Requests for test(s) where the RCI Laboratory acts as a Hospital Blood
	Bank must be made by a registered medical practitioner or an
	appropriately qualified healthcare professional acting on the
	instructions of a medical practitioner.
3.	Requests for referral services by Hospital Blood Transfusion
	Laboratories may be made by a medical scientist.
4.	Request for tests not processed in the RCI Laboratory are referred to
	specialist external laboratories. See Section 7 and section 7.2 for further
	details.

4.2 RCI Laboratory Request Forms

The laboratory service request must be accompanied by duly completed & <u>legible</u> RCI Laboratory Request form. A number of different request forms are available. These are used as outlined below.

Referral Service for Hospital Blood Transfusion Laboratories

BT - 0345 Request for Red Cell Immunohaematology Investigation Form

This form is used by referring Hospital Blood Transfusion Laboratories when they require the RCI Laboratory to undertake specialised investigations e.g. to resolve blood group serological anomalies, perform antibody investigations; or to request compatibility testing on these patients. This form may also be used by hospital blood transfusion laboratories to request ante-natal investigations, e.g. antibody titre and Anti-D/c quantitation.

Hospital Blood Bank Service Forms

BT - 0007 Blood Group and Compatibility Request Form

This form is used solely by facilities to which the laboratory provides Hospital Blood Bank services.

IBTS/RCI/FORM/0001 Request for Transfusion Reaction Investigation

This form is used by organisations, to whom the IBTS laboratory, acting as a Hospital Blood Bank has supplied compatible blood, when they wish to report a suspected transfusion reaction and request investigation of same.

BT- 0597 Haemovigilance Clinical Review Form

This form is used by the Haemovigilance Officer in facilities where the RCI laboratory acts as a Hospital Blood Bank, to document the patient information obtained following an adverse reaction or event. It should be forwarded to the IBTS SpMO/Registrar/Consultant Haematologist. It should include the details of the reaction/event, other relevant clinical information and results of haematology, biochemistry and microbiology tests performed as part of the adverse reaction investigation.

4.3 Ordering IBTS Request Forms

Request Form BT - 0345 **'Request for Red Cell Immunohaematology** Investigation' can be printed from the giveblood.ie website (<u>https://www.giveblood.ie/Clinical-Services/Red-Cell-Immunohaematology-Diagnostics/RCI-Test-Request-Forms/</u>)

All of the other forms are available on request from the RCI Laboratory by contacting personnel in the department.

4.4 Completion of Request Forms

A request form must accompany all samples referred for testing. Adequate completion of requests should include clinical information (e.g. obstetric history, transfusion history, reason for transfusion) so that work may be prioritised and processed accordingly in the laboratory; and to facilitate accurate result interpretation. As per BSH Guidelines the following mandatory patient personal identifiers must be provided on the request form and must be documented in a legible manner to be accepted for testing:

- 1. Patient's Surname
- 2. Patient's Forename
- 3. Patient's Date of Birth
- 4. Hospital number *
 - * Where the patient does not have a hospital number e.g. sample being referred from a GP (and the sample is not for compatibility testing purposes); an address will suffice as a third patient identifier in place of the hospital number.

The following information should also be documented on the request form:

- 5. Patient's gender
- 6. Patient's ethnicity
- 7. Location [referring hospital and ward (if given)]
- 8. Patient's address (* mandatory requirement if hospital number not applicable see above)
- 9. Details of the requesting clinician (& their contact details)
- 10. Date and time of sample collection (This is required for Hospital Blood Bank requests only)
- 11. Test(s) required
- 12. Number of units of blood required and date/time required (if for crossmatching)
- 13. Specific transfusion requirements for individual patients i.e. requirement for CMV negative and / or Irradiated blood

Where possible please provide the following information:

- 14. Relevant clinical information appropriate to the test(s) requested (e.g. clinical condition, medication)
- 15. Transfusion history (including results of serological investigations obtained by the referring centre, details of date of last transfusion, most recent haemoglobin level, historical antibodies, transplant history)
- 16. Antenatal history (including details of expected delivery date, anti-D administration, history of haemolytic disease of the foetus and new-born, history of intrauterine transfusions)
- 17. The specific clinical indication for a transfusion request
- 18. A clear indication as to whether the tests/services requested are urgent or routine

The Declaration(s) Must Be Signed:

- 19. The declaration regarding the correct labelling of the sample/request form and its validity <u>must</u> be completed (signed) by:
 - the person who took the sample (when the laboratory is acting as the sites Hospital Blood Bank)
 - the person referring the sample (Hospital Blood Transfusion Laboratory Referrals) Failure to complete the declaration may result in the sample not being processed.

Note:

Requests must be telephoned in advance if the service requested is urgent.

5 SAMPLE INFORMATION

Je.



sfective 2th Where the RCI laboratory acts as Hospital Blood Bank;

- For group and antibody screening one sample will suffice.
- For crossmatch requests, two samples collected at different times are • required, where there is no known historical ABO group, unless a secure bedside electronic patient identification system is in place. This is a recommendation from the BSH Guidelines 2012 to avoid transfusion of ABO incompatible blood due to misidentification of the patient at the time of sample collection. The exception is an emergency requirement for blood where one sample will suffice so as to not to unduly delay the transfusion. In this case group O blood will be selected for non-group O patients. It is important that the second sample is taken prior to transfusion so that the ABO/RhD group can be confirmed on the patient's cells without interference by transfused cells.

Where the RCI laboratory provides a Referral Service to Hospital Blood **Transfusion Laboratories**

Where RCI laboratory is providing a referral service, one sample is sufficient; however the responsibility for checking the historical group will reside with the referring Hospital Blood Transfusion Laboratory. If no historical group is available then the referring Hospital Blood Transfusion Laboratory should ensure the patient's ABO/RhD group has been verified on two separate samples prior to blood product issue.

Note: Re referral of compatibility requests to the RCI Laboratory

On a case-by-case basis and following discussion with the RCI Laboratory, the referring Hospital Blood Transfusion Laboratory may send segment(s) from suitable unit(s) for compatibility testing in conjunction with a patient sample(s). This is to expedite the provision of blood to the patient or where units of a particular phenotype are required and are already available from the referring hospital's blood

stocks. N.B. Segments must be labelled with the ISBT no. of the donor unit.

5.2 Service Requestor Responsibilities

- Obtaining consent from the patient for the tests required at the RCI Laboratory.
- Positively identifying the patient from whom the sample is taken.
- Safely disposing of the materials used in the collection of samples.
- Ensuring that samples containers meet the labelling requirements of the RCI Laboratory and that the request form has been completed to an acceptable standard.
- Ensuring that the test / services requested are appropriate.
- Ensuring that samples are delivered to the RCI Laboratory within a timeframe appropriate to the nature of the tests requested.
- Ensuring that appropriate transport containers are used (for the safety of all handlers).
- Ensuring that patient confidentiality is maintained.
- Ensuring that if referring unit segment(s) that the segment(s) are labelled with the unit ISBT No. and that the units meet the requirements of the patient's transfusion protocol.

5.3 Sample Labelling

5.3.1 Mandatory Requirements

The following essential information is MANDATORY on all samples referred to the RCI Laboratory and should be documented in a legible manner on the sample container:

- 1. Patient's Surname
- 2. Patient's Forename (initials are not acceptable)
- 3. Date of birth
- 4. Hospital number *
 - Where the patient does not have a hospital number e.g. outpatient / antenatal GP referrals (and the sample is not for compatibility testing purposes); an address/partial address will suffice as a third patient identifier in place of the hospital number
- 5. Date (and time where blood is requested [¥]) of sample collection
- 6. The initials/ signature of the person collecting the sample or in the case of an electronic labelling system, the label must contain unique identifier that can be used to trace the sample taker.
- ¥ Where the time of sample collection is not provided the sample time will be registered at the RCI Laboratory as 00:00 on the date of collection indicated.

Note:

All patient samples and forms must be labelled at the bedside applying positive patient identification.

5.3.2 Labels on Sample Tubes

- Sample tubes must never be pre-printed or pre-labelled.
- The Service Requestor's responsibility is to ensure that all printed labels for samples for blood transfusion testing are generated at the bedside and are compliant with BSH Guidelines (The administration of blood components: a British Society for Haematology Guideline, 2018).
- **Only labels that are printed 'On Demand'** next to the patient and immediately attached to the sample tube at the time of phlebotomy by the individual who took the sample are acceptable.
- Labels pre-printed away from the bedside or taken from the patient's notes (e.g. addressograph labels) are not acceptable on samples for processing. A repeat sample will be required.

5.3.3 Sample / Request Form Acceptance / Rejection

RCI Laboratory staff follow written standard operating procedures for the receipt and incoming inspection of samples and request forms. This is to ensure that samples taken for laboratory analysis can be accurately and unambiguously identified and that all necessary information is supplied for appropriate and timely analysis, interpretation and reporting.

Where the requirements with respect to labelling of the request form/sample container or sample quality issues are not met, this may result in the rejection of the request or a delay in sample processing.

Samples are accepted for testing if they are:

- 1. Of appropriate sample type for the tests required
- 2. Of sufficient volume for testing
- 3. If the information on the request form and sample are correctly matched
- 4. The sample & request card meet the mandatory labelling requirements.

Samples may be rejected in the following circumstances:

- 1. They are of an inappropriate sample type
- 2. They have leaked in transit
- 3. They are insufficient for testing
- 4. They are grossly haemolysed
- 5. They have been separated prior to referral
- 6. The sample and request form are mismatched, or the information is not correct
- 7. There is insufficient information on the sample and/or the request form.

- There is significant delay in receipt of sample from date/time of collection resulting in sample invalidity/instability. Note:
 Samples must be ≤ 7 days old on receipt with the exception of samples which fall under 72hr rule.
- 5.3.4 Non-Conforming Samples / Request Forms or Sample Quality Issues If a sample/request is identified as unacceptable, the referring laboratory/location or requestor (as appropriate) will be contacted and advised of any required corrective action or the need for a repeat sample in accordance with laboratory SOP. On occasion, rejected samples may be tested (see Section 5.3.5 Exceptions). In these instances, results reported will bear an appropriate caveat indicating the nature of the problem. A report will be generated for all rejected samples stating the reason for rejection. Sample rejection figures are monitored in the laboratory, if repeated issues are identified the referring hospital will be notified.

5.3.5 Exceptions

Exceptions may be made for samples from the following groups:

- Trauma, unconscious, or Emergency Department patients where the identity is not yet established. It is the responsibility of the referring laboratory to have a procedure in place for labelling of samples of unidentified patients. Ideally the minimum clinical information supplied should include: (1) a unique number, (2) gender and (3) approximate age. It is helpful to be informed of the ethnicity of the patient. Samples will be registered using the patient details on the sample tube. The sample details on the sample tube and request form must match.
- Where a repeat sample would be difficult to obtain and the result of testing is not to be used for transfusion purposes.
- Where the delay in acquiring a new sample might seriously prejudice a successful clinical outcome.
- Where the sample cannot be replaced, e.g. pre transfusion samples post transfusion reaction, samples taken at specific time periods e.g. foetal samples.

In the above exceptional circumstances, non-compliant samples may be accepted for testing with a documented authorised concession (e.g. written confirmation from the requestor verifying the patient identity) where delay in acquiring a new sample may seriously prejudice a successful clinical outcome for a patient, or where the sample cannot be replaced. In such cases the IBTS will not be responsible for errors made as a result of unacceptable labelling and/or samples issued by the referring facility. This may impact on the labelling and release of the suitable component, such that the



IBTS/RCI/CM/0001	Ver. 4	Page 22 of 51	
------------------	--------	---------------	--

component will be issued for transfusion at the discretion of the patient's clinician.

The decision to process the sample may require approval by the consultant at the IBTS. In all those instances the test report will identify and reflect the non-conforming issue.

5.3.6 Timing of Sample Collection Samples for compatibility testing should be referred to the laboratory without delay to facilitate timely testing of the samples and processing of requests.

Transfusions or pregnancy may stimulate the production of unexpected antibodies through either a primary or secondary response. The timing of samples selected for crossmatching or antibody screening must take account of this.

	Patient Category	Sample to be taken not more than
	Patient transfused or pregnant in the last 3 months	72 hours before transfusion ¹
	Patients receiving Anti-CD38 (with current or historical allo- antibodies)	
70	Patient not transfused or pregnant in the last 3 months Patients receiving Anti-CD38 therapy (no current/historical allo-antibodies)	7 days before transfusion ¹
	On-going cases	A formal deviation from the 3 day rule may be considered for patients that are being repeatedly transfused (e.g. AIHA, Myelodysplastic Syndromes) and have not become allo-immunised (i.e. have not formed clinically significant alloantibodies) allowing samples to remain acceptable for up to 7 days. This decision is agreed between the RCI Consultant Haematologist & the referring Hospital Consultant Haematologist, where a transfusion management plan is agreed.

Guidelines for the Collection of Samples

¹ This is the time between the sample being taken and the subsequent transfusion

5.3.7 Sample Storage

Whole-blood samples will deteriorate over a period of time. Problems associated with prolonged storage include red cell lysis, bacterial contamination, loss of complement in serum and decrease in potency of red cell antibodies, particularly IgM class antibodies.

5.3.8 Guidelines for the Storage of Samples: Pre-testing

BSH 2012 recommended working limits for the storage of blood testing samples (pre-analysis) are detailed below:

Patient Category	18−25 º C	2-8º C	-30º C
Patients transfused or pregnant in	Up to 48 hrs	Up to 3 days ¹	NA
the last 3 months	CC N		
Patients not transfused and not	Up to 48 hrs	Up to 7 days	3 months
pregnant in the last 3 months			

¹ This is the time between the sample being taken and the subsequent transfusion

5.3.9 Routine Referrals

2

Where samples are not being referred to the laboratory on the date collected they should be refrigerated at 2-8°C prior to transport unless otherwise specified in Section 7.1.1.

6 SAMPLE DELIVERY, PACKAGING AND TRANSPORT

6.1 Sample Delivery

Samples will be accepted by RCI laboratory at any time. They should be delivered to:

• Security at the National Blood Centre

Refer to Section 3.1 for map & picture of locations.

Note:

THE RCI LABORATORY MUST BE TELEPHONED IN ADVANCE OF URGENT REQUESTS AND THE SAMPLES DELIVERED AS PROMPTLY AS POSSIBLE.

6.2 Sample Packaging and Transport

It is advised that universal precautions be taken in the collection, packaging and delivery of the sample to the IBTS and that the patient's confidentiality is protected.

6.2.1 International carriage of dangerous goods by road

It is the responsibility of the service requestor to ensure the packaging; labelling and transportation of all samples comply with current European Agreement concerning Carriage of Dangerous Goods by Road Regulations. Legislation requirements are available

IBTS/RCI/CM/0001	Ver. 4	Page 24 of 51
------------------	--------	---------------

from the Health & Safety Authority website www.hse.ie. The requirements stated below apply to all diagnostic samples directed to the RCI laboratory.

6.2.2 Universal Packaging Procedure for the Transport of Diagnostic Samples

- Samples to be sent should be stored in a secure (preferably plastic) primary container.
- Wrap the sample tube/container in tissue or cotton wool which will act as absorbent material in the event of spillage.
- Place the sample tube/container in a biohazard bag.
- Place the biohazard bag with the sample tube and the request form in a padded envelope or an approved transport container.
- Label the envelope with a hazard warning label, "Diagnostic Sample, Category B UN 3373".
- Place the name, address and contact number of the destination laboratory on the outside envelope. *Note: It is very important to ensure that the address is correct and complete to ensure delivery to the correct location.*
- Address labels for samples boxes being referred to the RCI Laboratory are available at: - <u>https://www.giveblood.ie/Clinical-</u> <u>Services/Red-Cell-Immunohaematology-Diagnostics</u>
- Samples should be forwarded to the laboratory as soon as possible to preserve the integrity of the sample.
- Where blood is required the same day or it is an URGENT request, samples must be sent directly to the laboratory (see Note 2 below).
- The sample and the request form should be packaged so as to ensure patient confidentially at all times during transportation.

Note 1:

There is no requirement for a licensed courier to transport non-infectious diagnostic samples; however, An Post prohibits the sending of diagnostic samples by regular post.

Note 2:

Please contact the laboratory regarding all urgent samples. Ensure the transport box for urgent samples is marked 'Urgent'.

IBTS/RCI/CM/0001	Ver. 4	Page 25 of 51

7 TESTING PROVIDED

The table below lists the test procedures provided by the RCI Laboratory, tests that are further referred and other services available to customers. For further information on referral centres and consultants, refer to section 7.2.

Test / Service
RCI Laboratory
Antenatal Antibody Titration
Antibody Investigation
Investigation of Autoimmune Haemolytic Anaemia (AIHA)
Blood Group/Antibody Screen
Blood Group/Compatibility Testing
Blood Group/Compatibility Testing for Patients with Red Cell Antibodies
ABO Blood Group Anomaly Investigation
RhD Blood Group Anomaly Investigation (Serological)
Direct Antiglobulin Test
Elution
Extended RBC Phenotyping
Investigate Monoclonal Antibody Interference
Investigation of Haemolytic Disease of the Foetus and New-born (Where maternal red cell
antibodies are implicated / suspected)
Transfusion Reaction Investigation (Laboratory acting as Hospital Blood Bank)
Transfusion Reaction Investigation (Referred Sample)
Anti-D Quantitation
Anti-c Quantitation
Investigation of red cell polyagglutination
N

Samples may be further referred to other centres listed below

IBGRL: for complex Immunohaematological Investigation

NHSBT: for Investigation of IgA Deficiency & IgA Antibodies

NHSBT: for Cold Agglutin Titre and Thermal Range

NHSBT: for Paroxysmal cold hemoglobinuria (PCH) / Donath – Landsteiner

Other Services

Clinical & Scientific Consultancy Services

Haemovigilance Clinical Advisory Services (Hospital Blood Bank Service ONLY) Provision of HBB services including traceability and haemovigilance

IBTS/RCI/CM/0001	Ver. 4	Page 26 of 51
------------------	--------	---------------

7.1 Services Provided

- Pre-transfusion Compatibility / Specialised Immunohaematological Testing
- Provision of Rare Donor Red Cell Components from International Rare Blood Programmes.
- Concessionary Release of Blood Components
- Medical and Scientific Consultancy Service
- Haemovigilance Advisory Service
- Samples may be further referred to other centres where the testing is not performed in RCI or additional testing is required to support findings.

7.1.1 Pre-Transfusion Compatibility and Specialised Testing

The table that follows provides details of the tests available at the RCI Laboratory, sample requirements, any special requirements and turnaround times for test results.

• Initial verbal reports will be provided where emergency testing is performed or critical results are being reported.

Note 1

Turnaround time is defined as the time from sample reception at the IBTS to the time results/products are available for issue.

Note 2

Jerify when

Tests marked with an '*'are available out of hours for clinically urgent orders.

e 2ª				
IBTS/RCI/CM/0001	Ver. 4	Page 27 of 51		

Veritywhenin	Test Profile	Sample type (fresh venous sample)	Recommended Sample volume	Service details and requirements	Turnaround time test
	Antibody Titration	EDTA Whole Blood (WB)	1 x 6 ml	Batched testing For urgent testing contact the lab in advance	5 working days
	Antibody investigation * Including anti-G investigation	EDTA (WB)	2 x 6 ml	Progress will be discussed with requestor by telephone	5 working days
	Investigation of Autoimmune Haemolytic Anaemia *	EDTA (WB)	2 x 6 ml	Progress will be discussed with requestor by telephone	5 working days
	Blood Group / Antibody Screen ROUTINE or EMERGENCY *	EDTA (WB)	3 – 6 ml	Routine requests processed on next scheduled batch (See section 3.5)	Routine: Results usually available in 24 hrs.
	(Where RCI act as a Hospital Blood Bank)			Emergency: Processed immediately on receipt	Report will follow in 5 working days Emergency: ASAP - within 2 hrs of receipt of sample if

	2 ^A Y						
	IBTS/RCI/CM/0001		Ver. 4		Page 28 of 51		
Test Profile Sample type (fresh venous sample)		type (fresh venous	Recommended Sample volume	Service details and requirements		Turnaround time test	
				Contac labora advanc	tory in	no antibodies detectable	
	Blood Group and Compatibility Testing ROUTINE or EMERGENCY *	EDTA (WB)	1 x 6 ml	process batch Urgent requess Process immed receipt the IBT labora advand Segme suitabl be sen sample expedi provisi or whe particu	ised diately on t Contact TS tory in ce ents from le units may it with the e to ite the ion of blood ere units of a	Routine: 2-6 hours Urgent/emergency: ASAP (Within 2 hours of receipt of sample)	

Jerity when in

e 2 A Y				
IBTS/RCI/CM/0001	Ver. 4	Page 29 of 51		

Test Profile	Sample type (fresh venous sample)	Recommended Sample volume	Service details and requirements	Turnaround time test
			required & already available from the referring hospitals blood stocks N.B. Segments must be labelled with the ISBT no. of the donor unit.	
Blood Group and Compatibility Testing (Patients with red cell antibodies) ROUTINE or EMERGENCY *	EDTA (WB)	2 x 6 ml Minimum	Contact the RCI laboratory in advance Progress can be discussed by telephoning the RCI laboratory Refer to Section 3.5.2 for additional information. (See above re unit segments)	2-6 hours <u>Please note</u> this is dependent on the complexity of antibodies detected.
ABO Blood Group Anomaly	EDTA (WB)	1 x 6 ml	Telephone in advance if blood is	5 working days

Verifywhenin

e 2 ^{de re}				
IBTS/RCI/CM/0001	Ver. 4	Page 30 of 51		

1	Test Profile	Sample type (fresh venous sample)	Recommended Sample volume	Service details and requirements	Turnaround time test
	Investigation (Serological)			required for patient	
Verifywhenin	· · · · · · · · · · · · · · · · · · ·	EDTA (WB)	3 ml Cord Blood Sample	Serological testing will be performed on cord/neonate (<72 hrs old) samples only to determine maternal requirement for RhD Prophylaxis Progress will be discussed with requestor by telephone All other requests for RhD anomaly investigation should be referred to the Blood Group Genetics Laboratory at the NBC	Results phoned within 1 working day

		Q	24.40		
	IBTS/RCI/CM/0001		Ver. 4	Pag	ge 31 of 51
	Test Profile	Sample type (fresh venous sample)	Recommended Sample volume	Service details and requirements	Turnaround time test
Jerify when in	Direct Antiglobulin Test *	EDTA (WB)	1 x 6 ml	Next scheduled batch (See section 4.6) Note: there is no commercial IQC available for IgA and IgM portion of the monospecific DAT.	5 working days Where antibodies are investigated please see above
	Elution *	EDTA (WB)	1 x 6 ml	An eluate is only warranted if the patient has been transfused within the last month or there is evidence of haemolysis (or a delayed haemolytic transfusion reaction). Telephone in advance if blood is	5 working days

	IBTS/RCI	IBTS/RCI/CM/0001			Page 32 of 51	
	Test Profile	Sample type (fresh venous sample)	Recommended Sample volume	Service det and requir		Turnaround time test
				required for patient	or	
verify when in	Extended RBC Phenotyping	EDTA (WB)	1 x 6 ml	Extended phenotying recommen transfusion dependant patients an patients wi complex re antibodies. suitable for serological phenotypir patient mu have been transfused the previou months.	ded for n id ith ed cell . To be r ng the ist not within	5 working days
	Paternal Phenotyping	EDTA (WB)	1 x 6 ml			5 working days

IBTS/RCI/CM/0001	

Ver. 4

2440

0

,					
	Test Profile	Sample type (fresh venous sample)	Recommended Sample volume	Service details and requirements	Turnaround time test
Verify when in	Investigate Monoclonal Antibody Interference *	EDTA (WB)	2 x 6 ml	Progress will be discussed with requestor by telephone	5 working days
	Investigation of Haemolytic Disease of the New-born where maternal antibodies are implicated / suspected *	EDTA (WB) (Mother's sample) (Baby's sample)	1 x 6 ml 1-3 ml	Note: Investigations will be limited to Red Cell Serological studies Processed in next scheduled batch (See section 4.6) Must be telephoned in advance	Results phoned within 1 working day from availablility of result (Quantitation/ Antibody Titration)
	Transfusion Reaction Investigation * (Where the IBTS acts as a Hospital Blood Bank)	EDTA (WB) Post transfusion + Clot activated post	2 x 6 ml 1 x 6 ml	Must be telephoned in advance. Contact medical consultant / medical registrar on duty / on call, for direction Please return implicated unit (if	ASAP 2-6 hours of receipt of sample for initial serological results Note: Where bacteriological screening of the implicated units is

		.0	24			
	IBTS/RCI/	IBTS/RCI/CM/0001		Ра	Page 34 of 51	
	Test Profile	Sample type (fresh venous sample)	Recommended Sample volume	Service details and requirements	Turnaround time test	
Verifywhenin		transfusion sample The implicated unit must be sealed by a coupler and returned		available) and the administration set (if possible). (Even an 'empty pack' may provide a sample from an attached segment) The remaining un- transfused units must be quarantined at the hospital or returned to the IBTS, pending medical release. Part B (white) of the traceability label (BT396) must not be removed from the units when returning to the IBTS.	required, or immunological investigation is necessary, the turnaround time may be extended beyond 7 days A written report of the serological results only may be available within 5 working days	

IBTS/RCI/CM/0001	

Ver. 4

e 2ª ro

	Test Profile	Sample type (fresh venous sample)	Recommended Sample volume	Service details and requirements	Turnaround time test
Jerity when in	Transfusion Reaction Investigation * (Referred samples from Hospital Blood Transfusion Laboratories)	EDTA (WB) Pre Transfusion (if crossmatch performed by referring hospital) Post transfusion + A clotted post transfusion	1 x 6 ml 2 x 6 ml 1 x 6 ml (if required)	Must be telephoned in advance. Progress will be discussed with requestor If pack culture is to be performed by the IBTS, the implicated unit must be sealed by a coupler and returned. Note: Testing may also be performed	ASAP 2-6 hrs of receipt of sample for initial serological results. A written report of the serological results only may be available within 5 working days
	Anti- D Quantitation	sample should also be referred where possible EDTA (WB)	1 x 6 ml	on the post- transfusion sample only, as requested by referring hospital. Batched testing: Performed 2-3 times weekly	Result phoned within 5 working days


	2 A Yo							
	IBTS/RCI/C	Ver. 4	Ver. 4 Pag		e 36 of 51			
	Test Profile	Sample type (fresh venous sample)	Recommended Sample volume		e details equirements	Turnaround time test		
1.				Separated samples will not be processed		Written report 7 working days		
	Anti-c Quantitation	EDTA (WB)	1 x 6 ml	Batched testing: Performed weekly Separated samples will not be processed		Result phoned within 5 working days Written report 7 working days		
	Investigation of red cell polyagglutination	EDTA (WB)	1 x 6 ml	give ris groupi anoma associa necrot enterco paedia These usually transfi	ination may se to blood ing alies and is ated with	5 working days		

Verify when in

2AXC						
	IBTS/RCI/CM/0001		Ver. 4		Page 37 of 51	
	Test Profile	Sample type (fresh venous sample)	Recommended Sample volume		e details equirements	Turnaround time test
Jerity when in					onent on is icated t resolution erred sample	
7er.				advan	atory in	

* Testing performed during both routine and out of hours service

Note: With regard to sample volume and number of samples required, exceptions may be made for patients where blood is difficult to obtain such as neonatal or paediatric patients.



7.2 Referral Test Services

The IBTS act as the national facilitator for the referral of samples to NHSBT Barnsley for the investigation of IgA deficiency, IgA antibodies and Cold Agglutinins. In addition following investigation by the IBTS, samples may be requested for referral to the IBGRL for confirmatory testing or where the results obtained by the RCI laboratory are inconclusive. The NHSBT user guide may be accessed at http://hospital.blood.co.uk./diagnostic-services/user-guides/. The IBGRL user guide may be accessed at <a href="http://ht

All referral laboratories listed below are accredited to ISO 15189 for the tests described.

Note 1

Referrals to the external laboratories attract testing and transportation charges. Additional time will be incurred where samples are referred externally. In such cases the scientific staff at the IBTS will inform the requestor.

Note 2

The report from the external laboratory will be forwarded to the referring hospital (a copy of this report will be retained by the IBTS for reference).

IBTS/RCI/CM/0001 Ver. 4 Page 39 of 51

 $\langle \rangle$

Test profile /service	Centre	Sample type (fresh venous sample)	Sample volume	Service details and requirements	Turnaround time test
Complex Immunohaematology Testing	IBGRL	EDTA (WB)	Contact RCI laboratory prior to sending.	Samples are accepted from overseas reference laboratories ONLY. Service includes confirmation of rare specificities previously determined by the RCI laboratory and determination of possible underlying specificities. Larger samples are preferable and sometimes may be essential. Anti-coagulated samples should not be separated. Urgent referrals are defined as those where blood for transfusion is needed as quickly as possible. The RCI laboratory will contact the IBGRL by telephone to discuss reason for referral.	Due to the varying nature of the requests sent to Red Cell Reference a turnaround time cannot be specified. The time between receipt of sample and reporting will depend on the clinical situation of the patient, the complexity of the investigation and the number of samples in the laboratory at any given time. Cases are prioritised and therefore some investigations may be necessarily delayed.

	IBT	S/RCI/CM/0001	çi ⁿ	Ver. 4	Page 40 o	f 51	
Investigation of IgA Deficiency & IgA Antibodies	NHSBT Barnsley	EDTA (WB)	Contact RCI laboratory prior to sending.	indications Samples for in	vestigation should the RCI laboratory eferral to the	receipt Repeat testi extension to	ays from sample ng will incur an the turnaround ssibly may require ling
Jerify when in							

			24				
	IBTS/RCI/CM/0001 Ver. 4 Page 41 of			f 51			
Cold Agglutinins/CHAD Investigation (Investigation comprises DAT, room temperature screen, cold titre, thermal amplitude as necessary.)	NHSBT Barnsley	EDTA (WB) Serum sample * (separated at 37C)	Contact RCI Laboratory prior to referring samples.	A DAT is perfor plasma is inve- presence of cl red cell alloan Subsequent to performed at indicated to e significance at amplitude. *It is not nece warm separate titration studi specifically re the RCI labora primary samp separated sam labelled with identifiers	prmed and the estigated for the linically significant atibodies at 20C. ests may need to be 4C and 30C where stablish clinical nd thermal essary to re samples unless es are required or quested to do so by atory. Send the le tube of the nple tube that is the patient	Result availa days of samp NHSBT. Report will b	ble within 5 working ble receipt by e despatched by the bry following receipt
				Cold	med on request in in Disease patients.		

	$e^{2^{A}}$							
	IBT	S/RCI/CM/0001			Ver. 4	Page 42 of	51	
Paroxysmal cold	NHSBT	EDTA	Conta	act RCI	Biphasic haem	olysins as a cause	Result availa	ble within 5 working
hemoglobinuria (PCH)	Barnsley		Labo	ratory	of AIHA are ex	tremely rare and	days of samp	ole receipt by
/ Donath –		&	prior	to	mainly seen as	s a post-viral event	NHSBT.	
Landsteiner			refer	ring	in children. Ro	outine investigations		
		Separated	samp	les.	for AIHA do no	ot include the test	Report will b	e despatched by the
		serum sample			for biphasic		RCI Laborato	ory following receipt
		which has			haemolysins b	ut where indicated,	of same.	
		been			or on request,	the Donath-		
		separated from			Landsteiner te	est can be		
		a whole blood			performed if p	aroxysmal cold		
		sample that			haemoglobinu	ıria (PCH) is		
		has been			suspected. If p	oositive, the		
		allowed to clot			specificity of t	he antibody can be		
		at 37oC.			determined to	confirm the		
					diagnosis.			

7.3 Provision of Rare Donor Red Cell Components from International Rare Blood Programmes

Where the patient requires red cells of a specific red cell antigen profile that is not available in the IBTS stock supply or on the IBTS donor panels, where appropriate a request will be made to an International Rare Blood Bank Programme for the required number of units (either from current stock, following donor call up or frozen blood stocks).

Each individual case will be discussed with the IBTS medical staff and authorised on a Consultant to Consultant basis, to determine the exact requirements for individual patients and advise of associated difference in transfusion risk profile as appropriate. Procurement of product is dependent on the availability of the blood. Walk in whole blood donation and liquid red cell component transfusion is preferred and this requires scheduling. In clinical emergencies frozen recovered red cell components may be sourced Approval of medical consultant is required. Please provide the maximum notification possible for this service.

Note:

Direct consultation with the medical consultant on duty will be required.

7.4 Concessionary Release of Blood Components

Concessionary release of blood components, or acting contrary to an SOP, is sometimes the necessary and appropriate course of action in the best interest of patients. To act contrary to an SOP requires prior authorisation, or justifiable authorisation as soon after as is practicable, by the IBTS Consultant Haematologist or other suitably competent person who should discuss the clinical consequences with the clinicians in charge of the patient this process is proceduralised in IBTS/MED/SOP/0050.

Conditions which require concessionary release procedure:

- Use of D positive blood for a D negative patient who would normally be excluded from receiving D positive units.
- Use of antigen positive or un-typed red cells in patients with atypical red cell antibodies.
- Issue of red cells to patients with autoimmune haemolytic anaemia (AIHA) without the necessary exclusion of underlying antibodies.
- Issue of components that do not meet known special requirement, e.g. CMV negative or irradiated.
- Where it is necessary to act contrary to a Standard Operating Procedure in the best interest of a patient, this will be handled in accordance with relevant IBTS Quality Assurance documentation.

The event will be recorded on a Concessionary Release form. The name and designation of the requesting clinician who has agreed to accept the concession for the patient will be recorded along with the details of the IBTS Clinician who

IBTS/RCI/CM/0001	Ver. 4	Page 44 of 51
------------------	--------	---------------

has authorised the release of the blood product/component. A copy will be sent/faxed to the hospital blood transfusion laboratory; or to the requesting clinician where the IBTS laboratory provides Hospital Blood Bank services.

7.5 Medical and Scientific Consultancy Service

The IBTS will provide medical and scientific advice for all the above services. These services are available at all times with respect to blood transfusion practice. For contact names and numbers see Section 3.4 of this manual.

7.6 Haemovigilance Advisory Service

All haemovigilance queries relating to situations where the RCI laboratory acts as a facilities Hospital Blood Bank should be directed to IBTS medical staff, directly to laboratory senior scientific staff or to the IBTS Biovigilance Officer.

7.7 Repeat Examination

It is the policy of the laboratory, in the event of analytical failure to:

- Repeat the test using the relevant procedure(s) or
- Store the sample in appropriate conditions, until the cause of the analytical failure is identified and corrected; and then repeat the test. The urgency of the outstanding sample is reviewed by the relevant laboratory director or nominee.
- Samples are retained in accordance with local Guidelines for Storage of Examined Samples for Archive and Look Back Purposes.
- Should additional samples be required the laboratory will contact the requesting location.

7.8 Further Examination of the Primary Sample

Where further testing is relevant to the investigation, then it is the policy of the laboratory to pursue further investigation using the primary sample. If additional investigations / blood products are required please contact the laboratory to ensure that sufficient sample is available and that the sample is still valid. Red cell samples are stored at 4°C for seven days or until the report is issued. Separated serum/plasma samples at stored at -20°C or lower for 14 days. Antenatal samples; Anti-D, Anti-c quantitation and titre samples are held for parallel testing for up to 10 months.

8 EXTERNAL AND INTERNAL QUALITY ASSURANCE SCHEMES

8.1 External Quality Assessment Programmes (EQA)

The RCI Laboratory participates in relevant available external third party assessment schemes.

This includes schemes operated by:

- UK NEQAS (United Kingdom National External Quality Assurance Scheme) for Blood Transfusion Laboratory Practice
- UK AQQAS (Antibody Quantitation Quality Assurance Scheme) for Anti-D/c quantification
- Aurevia for Direct Antiglobulin Testing
- WASPS (Welsh Assessment of Serological Proficiency Scheme) for compatibility testing

External proficiency testing is performed by all staff working in the RCI laboratory on an annual frequency.

The laboratories are committed to participating in other schemes as they become available and are required to ensure comprehensive assessment of the test repertoire.

8.2 Inter-Laboratory Comparison Scheme

The RCI laboratory also participates in an inter-laboratory comparison scheme for elution, adsorption, titration, neutralisation, and DTT treatment test methods as no formal EQA programmes are available for these particular test methods (with the exceptions of pilot exercises).

8.3 Internal Quality Assessment Programme

Internal controls are included in all tests: no tests can be accepted or reported unless control results are acceptable.

All test procedures are covered by Standard Operating Procedures and only trained and authorised staff may perform procedures. Staff competency is also assured before a staff member may perform a procedure. All procedures are regularly reviewed.

8.4 Authorisation

The RCI Laboratory is part of the IBTS quality management system, which is covered by GMP (Good Manufacturing Practice) and is inspected yearly by the HPRA (Health Products Regulatory Authority), the Blood Establishment Authorising Authority. The IBTS Blood Establishment is Authorised under BE Number **0002**.

8.5 Non-Conformance / Failure in an EQA Scheme

Non-conformances are managed by controlled procedures, with investigation, corrective and preventative actions and review of practices, taken as

IBTS/RCI/CM/0001	Ver. 4	Page 46 of 51	
------------------	--------	---------------	--

appropriate. The laboratory may be audited at any time provided that the IBTS Director of Quality and Compliance and the Laboratory Directors are notified in advance and that the time is agreed by all parties. Where the RCI Laboratory fails an external EQA scheme, all customers of the service will be notified.

9 REPORTING OF RESULTS

9.1 Approval of Test Results and Issuing Reports

- All test results are reviewed and approved by a medical scientist before release.
- Valid results of automated testing are entered electronically into the IBTS computer system. Results of manual testing are entered by two medical scientists (or approved system during on-call hours). All results are validated by the IBTS computer system.
- Where relevant clinical advice and interpretative comments will be included on the test report.
- Clinical advice and interpretative comments are based on recommendation from BSH guidelines as standard.
- Where blood has been crossmatched and issued for a patient the units will be tagged with an IBTS compatibility tag and a hardcopy report will be issued with the blood.
- Routinely hard copy reports are printed and posted to the requesting laboratory/location. One hardcopy will be sent for each sample/request.
 - All reports are checked and signed by the Chief Medical Scientist or other senior person in charge once testing is complete.
 - Where blood has been issued for a patient the accompanying hardcopy report issued with the blood will be signed by the medical scientist who issued the blood.
 - Relevant reports are also reviewed and signed by the Consultant Haematologist.

Compatibility results may be reported as compatible, least incompatible or suitable in accordance with BSH Guidelines and Daniels et al. The term 'Suitable' printed on the Compatibility Report, indicates the units are compatible/least incompatible for patients with auto-antibodies using adsorbed plasma, this term is recommended in the BSH Guidelines 2012. This term is also used to report that units were compatible following monoclonal antibody investigation.

- Where the RCI laboratory has crossmatched segments from suitable units which were provided by the referring hospital, the unit number and segment numbers will be listed on the report and identified as compatible/least incompatible/ suitable as appropriate.
- Reports are also available to hospitals on demand (where urgently required) and are issued with all emergency requests.
- The laboratory retains a copy of the report with the original request form. Where a preliminary or interim report is issued, a final report will follow.
- Only scientific or medical staff may issue reports to the hospital. The report must be issued to the medical / clinical personnel responsible for

the patient or to scientific staff at the referring laboratory.

- Cumulative patient reports / worksheets for each patient are stored together.
- Where no component has been issued and an antibody investigation has been carried out by the RCI Laboratory an antibody report will be issued within 5 working days (Refer to Section 7.1.1). Should the report be accompanied either by a clinical comment or a covering letter from the director of the laboratory this will incur a further delay in the reporting times.
- A written report will be issued within 14 working days from the receipt of the sample in all cases; except where samples have been referred to an external laboratory.

9.2 Issuing Reports on Critical Samples where Results are Delayed.

It is policy to immediately notify the referring hospital / team when there are indications that the results may be delayed. The laboratory will maintain a record of all such correspondence.

A verbal report will be given as progress of the test becomes available, if required.

Note:

The RCI laboratory will not release results of examinations performed directly to patients. If requested, the laboratory will advise that the best practice is to issue the result to their clinician who can then discuss the implications of the test results with them.

9.3 Reporting of Results by Fax / Encrypted Email

Where clinically requested by the referring laboratory or the hospital clinician, the IBTS will issue results by fax or encrypted email. Where reports are requested to be faxed / emailed they will be signed by a scientist. They may not have a Chief Medical Scientist or Consultant Haematologist signature and may be labelled as 'preliminary' if testing has not been fully completed or as 'interim' if the scientist is lone working. Telephone verification of the receipt of the faxed results is required. A hard copy of the report will follow in the post. Reports may be sent by email on request, only where a secure encrypted email process has been put in place with the requestor.

9.4 Telephoned Results

The RCI Laboratory provides telephoned results (to the patient's clinician / designated clinical personnel or to the referring laboratory) as clinically required.

Criteria for telephoning results includes the following:

- 1. Significant unexpected findings
- 2. When there is a significant delay in turnaround time
- 3. When blood / blood components are ready for issue/delivery
- 4. Relevant antenatal testing results.

5. When requested by the referring location

When requesting a verbal report, the patient's personal identifiers i.e. patient's name, date of birth and hospital number must be given to the RCI scientific staff.

The RCI staff will also require the details of the requestor i.e. their own name.

In accordance with laboratory procedures a record of all verbal reports is maintained by the laboratory. A hard copy of the report will follow all verbal reports.

9.5 Archiving of Patients Records

It is IBTS policy to store copies of original request forms and the reports issued by the IBTS laboratories for >30 years, either by electronic or paper record systems.

10 CUSTOMER SERVICE / SATISFACTION AND REVIEW

10.1 Service Level Agreements

Customer satisfaction is assessed through a yearly survey of users, feedback received at Hospital Transfusion Committee meetings and processing of complaints. Customer complaints are reviewed and discussed at the RCI SMS Meetings and at the RCI Annual Quality Review Meetings.

10.2 Customer Complaints / Compliments

The IBTS customer complaint form (IBTS/QA/SOP/0063 Attachment 6.2) for blood products, defects & testing / reporting of patient samples is available of the giveblood.ie website at the following link:

https://healthprofessionals.giveblood.ie/clinical-services/qualitycompliance/quality-documentation/quality-documentation.html. The form is

Completed by the user and emailed to <u>qualityassurance@ibts.ie</u>. A response will be received following investigation of the complaint. Refer to IBTS/RCI/LM/0001 Section 7.7

The RCI laboratory welcome user feedback, this form is available from the following section of the website: <u>https://healthprofessionals.giveblood.ie/clinical-services/transfusion-transplantation/red-cell-immunohaematology-diagnostics/</u> On completion of the patient, user and personnel feedback form, it is emailed to <u>rci@ibts.ie</u>. The laboratory will respond to the feedback once received.

10.3 Quality Management Review

Refer to Section 8.9 Management Review of IBTS/RCI/LM/0001

10.4 Customer Liaison

10.4.1 Refer to Section 8.6.1 IBTS/RCI/LM/0001 for details on Customer Satisfaction.

IBTS/RCI/CM/0001	Ver. 4	Page 49 of 51	

10.4.2 Samples may be requested from the RCI laboratory for validation purposes. Please email requests to <u>rci@ibts.ie</u>.

10.5 Hospital Transfusion Committees

Refer to Sections 11 of IBTS/RCI/LM/0001 for Hospital Transfusion Committees.

10.6 Continuous Improvement

Refer to Section 8.6 Continuous Improvement of IBTS/RCI/LM/0001

11 TRACEABILITY AND REPORTING OF SERIOUS ADVERSE REACTIONS (SARS) AND SERIOUS ADVERSE EVENTS (SAES)OTHER REQURIREMENT

11.1 Traceability

- **11.1.1** Refer to Section 11 of IBTS/RCI/LM/0001 for additional information.
- **11.1.2** A traceability label, BT 396 is attached to all blood / blood components with a cable tie when issued from the RCI laboratory.

Jerity Marin Use.

STOP SEE BACK OF	THIS TAG BEFORE TRANSFUSION	
Irish Blood Transfusion Serv Seirbhís Fuilaistriúcháin Eire	NBC 01 4322800 Fax 01 4322930	PRE ADMINISTRATION
		A REAL PROPERTY OF THE REAL PR
Donation Component:		STEP 1: Check the component has been prescribed Check any special requirements e.g. irradiated Check if concomitant drugs prescribed e.g. diuretic.
Signature 1:	Date Given:	STEP 2: Check and document baseline observations.
Signature 2	Time Given:	STEP 3: Check expiry date and time of component.
B 1 51 1 1 1		Check pack for leaks, discolouration or clumping.
	I place in patient's Medical Records	
Surname:	Forename:	ADMINISTRATION
DOB:	Gender	STEP 1: Ask the patient to tell you their Surname,
		Forename and Date of Birth. Be especially vigilant
Hospital:		with unconscious or compromised patients,
		refer to your local hospital policy.
Ward:		STEP 2: Check their Surname, Forename and Date of
Hospital No:	Cutable for Tree for the state	Birth and Patient Identity Number against their
isopital NO.	Suitable for Transfusion Until	wristband and the compatibility label.
	For internal use only refer	STEP 3: Check that the information on the compatibility
	to patient report for expiry	label matches the details on the blood component
Patient's Blood Group: Con	nponent: Comments:	i.e donation number, blood group.
		If there are any discrepancies – DO NOT PROCEED - contact your Hospital Transfusion Laboratory and HVO.
Special Requirements / Trans	sfusion Protocol:	If you suspect a transfusion reaction- STOP the
		transfusion immediately, seek medical advice,
Donation Number:		and contact the HVO and Transfusion Laboratory.
completed section below	een started, you must send the back to the Hospital Transfusion policy. This is a legal requirement	
Burname:	Forename:	Under SI # 547 of 2006 European Communities
lospital No:		(Human Blood and Blood Components Traceability
iospitar No.	Lab Sample No:	Requirements and Notifications of Serious
Oonation Number:		Adverse Reactions) Regulations 2006.
onation Number:	DOB:	
		IT IS A LEGAL REQUIREMENT
omponent:		that this section of the label be completed and
		returned to the Transfusion Laboratory
Date Given:	Time Given:	
I confirm that the above name	d patient received this blood component.	© Irish Blood Transfusion Service BT396-2
ign and Print Name	Hosp.	Institution Transfusion Service D1390-2
	Wd.	

11.2 Serious Adverse Reactions (SARs) and Serious Adverse Events (SAEs) 11.2.1 Refer to Section 11.2 of IBTS/RCI/LM/0001

12 REFERENCES

- 12.1 Milkins, C. et al. (2013), Guidelines for pre-transfusion compatibility procedures in blood transfusion laboratories. Transfusion Med, 23: 3-35. [IBTS/EXT/DOC/0025].
- 12.2 Blood Directive Directive 2002/98/EC ~ "Setting the standards of quality and safety for the collection, testing, processing, storage and distribution of human blood and blood products and amending Directive 2001/83/EC". [IBTS/EXT/DOC/0012]
- 12.3 EU Directive 2004/33/EC Annex IV titled "Storage, Transport and Distribution Conditions for Blood and Blood Products". [IBTS/EXT/DOC/0012]
- 12.4 SI 360 / 05 European Communities (Quality and Safety of Human Blood and Blood Products) Regulations 2005. This is the statutory instrument which adapts the EU Directives as defined above Into Irish law. [IBTS/EXT/DOC/0012]
- 12.5 Traceability SI 547/06 Compliance with Article 14 (Traceability) and Article 15 (Notification of Serious Adverse Reactions and Events) of EU Directive 2002/98/EC". [IBTS/EXT/DOC/0012]
- 12.6 Minimum Requirements for Blood Bank Compliance with Article 14 (Traceability) and Article 15 (Notification of Serious Adverse Reactions and Events) of EU Directive 2002/98/EC (AML-BB). [IBTS/EXT/DOC/0017].
- 12.7 Directive 2005/61/EC ~ "Implementing Directive 2002/98/EC of the European Parliament and of the Council as regards traceability requirements and notification of serious adverse reactions and events". [IBTS/EXT/DOC/0012]
- 12.8 Directive 2005/62/EC ~ "Implementing Directive 2002/98/EC of the European Parliament and of the Council as regards community standards and specification relating to a quality system for blood establishments". [IBTS/EXT/DOC/0012]
- 12.9 ISO 15189. Medical Laboratories Particular requirement for quality and competence, 2022. International Organisation for Standardisation.
 [IBTS/EXT/DOC/0033].
- 12.10 NHSBT user guide may be accessed at https://hospital.blood.co.uk/diagnosticservices/user-guides/
- 12.11 IBGRL user guide may be accessed at https://ibgrl.blood.co.uk/services/red-cellreference/

13 ATTACHMENTS

- **13.1** List of customers for which RCI lab provide referral testing services
- 13.2 List of customers for which RCI lab provide HBB services

Customers for which RCI lab provide referral testing services	
Cork University Hospital	
Bon Secours Hospital Cork	
Mercy Hospital	
General Hospital Bantry	
Mallow General Hospital	
General Hospital Cavan	
General Hospital Letterkenny	
St.Vincent's PUBLIC Hospital	
St.James Hospital	
Tallaght Hospital	
Beaumont Hospital	
Mater PUBLIC Misericordiae Hospital	
Mater PRIVATE Hospital	
Temple Street Hospital	
Crumlin Hospital	
Coombe Women's Hospital	
National Maternity Hospital	
Rotunda Hospital	
Bon Secours Hospital Dublin	
Blackrock Clinic	
St.Mary's Hospital Cappagh	
James Connolly Memorial Hospital	
St.Lukes Hospital, Rathgar	
Beacon Hospital	
Hermitage Medical Clinic	
U.C.H. Galway	
Portiuncula Hospital	
Bons Secours Hospital Galway	
The Galway Clinic	
County Hospital Naas	
St.Lukes General,Kilkenny	
Kerry General Hospital	
Bon Secours Tralee	
Lady of Lourdes, Drogheda	
University Hospital Limerick	
General Hospital Portlaoise	
Our Lady's Hospital Navan	
General Hospital Castlebar	
General Hospital Tullamore	
Co.Hospital Roscommon	
General Hospital Sligo	
Sth Tipperary Gen. Clonmel	
Waterford Regional Hospital	
Midland Regional, Mullingar	
General Hospital Wexford	

	IBTS/RCI/CM/0001	Ver. 4	Attachment 13.1	Page 1 of 1
--	------------------	--------	-----------------	-------------

Customers for which RCI lab provide HBB services

Our Lady's Hospice, Harold's Cross

Royal Victoria Eye & Ear

Verimmen muse. status cupretime the etime a feature of the status cupretime the etime and the status cupretime a status cupreti Blackrock Hospice

IBTS/RCI/CM/0001	Ver. 4	Attachment 13.2	Page 1 of 1